

HPCAT  
16 IDD


---

# User Manual



*High Pressure Collaborative Access Team  
Advanced Photon Source*

# 16 IDD User Manual



# Table of Contents

<b>ABOUT 16 ID-D</b>	<b>1</b>
BEAMLINE DESCRIPTION	1
MAJOR ID-D FEATURES	2
<b>TROUBLESHOOTING</b>	<b>3</b>
LOST BEAM INTENSITY	7
BLUEDIAMOND	6
<i>Sluggish response</i>	6
<i>No response</i>	6
IDL SCANSEE SOFTWARE IS NOT WORKING	7
LOSING BEAM INTENSITY	4
<b>BEAMLINE OPERATIONS</b>	<b>5</b>
LOGIN	5
HOW TO START EPICS INTERFACE	5
BLUEDIAMOND	9
IDL SCANSEE SOFTWARE	11
2D SCAN OF SAMPLE ABSORPTION	13
STEERING BEAM PROCEDURES	16
DCM ENERGY CALIBRATION USING FOIL	23
OVERPLOT SEVERAL SCANS IN BULEDIAMOND	31
<b>SOFTWARE</b>	<b>37</b>
INTRODUCTION	37
SOFTWARE AT CONTROL COMPUTER	37
SOFTWARE AT DATA ANALYSIS COMPUTER	37
SOFTWARE AT FOCUS CHECKING COMPUTER	38

## About 16 ID-D

### **Beamline description**

The The HPCAT 16ID-D beamline is dedicated to x-ray scattering and spectroscopy research of materials under high pressure, typically in diamond anvil cells (DAC). With beams of 2meV energy resolution at  $^{57}\text{Fe}$  energy of 14.413keV, nuclear resonant techniques are used for determining phonon density of state and Mössbauer effect for Fe-containing samples at high pressures. With beams of 1eV energy resolution, a high efficient instrument for x-ray inelastic scattering (x-ray Raman) provides information of charge dynamics and chemical bonding at high pressures. The 16ID-D station has also a setup of emission spectrometer, measuring the details of fluorescence signals for information of spin state and structures of valence electrons at high pressures.

#### **Beamline scientists:**

Paul Chow

Phone: (630) 252-4015    E-mail: pchow@ciw.edu

Yuming Xiao

Phone: (630) 252-4058    E-mail: yxiao@ciw.edu

## Major ID-D features

Feature	Description
Source	Dual undulators Type A, single or in tandem mode
Monochromator	Oxford Diamond (1 1 1) DCM
Energy Range	5-25keV
Beam Size and Focusing Optics	35 (V)x 50 (H) $\mu\text{m}$ with 1 meter long KB mirrors 10 (V)x 10 (H) $\mu\text{m}$ with 200mm long KB mirrors or combine 1 long KB mirror with 1 short KB mirror
Flux at Sample Position	NRS: $\sim 1 \times 10^9$ ph/s, 2meV energy resolution Others: $\sim 2 \times 10^{12}$ ph/s flux based on using 1 meter long KB mirrors
Established Techniques	Nuclear Resonant Inelastic X-ray Scattering (NRIXS) Nuclear Forward Scattering (NFS) Inelastic X-ray Scattering (IXS) X-ray Raman Scattering (XRS) Resonance/X-ray Emission Spectroscopy (R/XES) Partial Fluorescence Yield X-ray Absorption Spectroscopy (PFY-XAS) X-ray Fluorescence Spectroscopy (XFS)
Detectors	100K Pilatus detector, Avalanche Photo Diode (APD), AmpTek
Support Equipments	Cryostats for XES and NFS, online Ruby system

## Troubleshooting

### Data collection computer crash

If the data computer crashes or doesn't response during the experiment, please restart the computer and follow steps in "LOGIN" and "HOW TO START EPICS INTERFACE" section in Chapter 3.

### Bluediamond

#### *Sluggish response*

The real-time data display is basically relying on the data in/out control routine which uses EPICS channel access method; therefore this software becomes sometimes sluggish depending on the local network traffic status, which may vary station by station even within the same sector and at the same time. Many bugs have been reported by users for past cycles; and most of them were fixed in the latest version (2010 March).

#### *No response and Restart*

If the software is waiting for the input PV values updated as a scan started and if there is no update in the PV, it may wait for the input signal update forever. There is no feasible timeout routine to reset this error, according to the writer of software. This error is, however, often caused by the errors involved in EPICS Scan procedures, which can be caused by a mistake in the input PV values (e.g., incorrect scan range) and may need to acknowledge the error or reset it on the MEDM interface (not on the Bluediamond).

If you need restart Bluediamond, firstly you need press "Ctrl+Alt+Del" to start the Task manager under Windows, then kill the Bluediamond in the Application. Secondly, follow steps in "BLUEDIAMOND" section in chapter 3.

## **IDL scanSee Software is not working**

If ScanSee crashes or has no response, please close it and restart it following steps in “IDL SCANSEE SOFTWARE” section in Chapter 3

## **Losing beam intensity**

One of the most frequent causes of beam intensity decrease in the experimental hutch is accidental bump into the pinhole, if so, please check the Y and Z positions of pinhole and the sample alignment.

If it is not the pinhole problem, there are several possibilities.

- 1) If the APS is running at decay mode (ring current decreases gradually from 100mA to 80mA), you will see intensity decrease slowly , that is normal, so you don't do anything. The intensity will come back after re-injection.
- 2) If there is a beam lost, the DCM need time to warm up, during warming up period, the intensity will drop, to get intensity back, we need adjust DCM pitch to get maximum intensity. After adjust pitch, beam position will also be changed, so we need steer beam positions, please follow steps in “STEERING BEAM PROCEDURES” in chapter 3.
- 3) If the intensity changes suddenly, you can look at APS beam position monitor to check source beam positions and compare it to those values we print out. If you see big change, please let beamline scientist know, then we need call floor coordinator to ask APS control room to steer beam.

## Beamline operations

### Login

Both beamline control computer and data analysis computer (Windows XP Professional operating systems) has identical login id and password for users:

login: **16idd\_user**

password: **idduser@16**

This login id and password will be required when re-logging in the computers (e.g., after rebooting the operating system) and also valid to access network data folder <\\hpcat21\16idd>.

### How to start EPICS interface




1. double click "MEDM"







2. change medm.exe from “Edit” to “execute”



3. click “IDA” “IDC” “IDD” from 16id\_main.adl



4 click “Monitors” from 16idd.adl



then scaler window pops up

scaler16\_more.adl


Done OneShot time Count time Elapsed time  
Count AutoCount 1.00 1.000 1.000

#	Description	Gate?	Preset count	Actual count
1	Timer	N Y	10000000	10000000
2	IDC ICO V1 -I	N Y	0	0
3	IDC ICO V2 -I	N Y	0	0
4	IDC ICO H1 -I	N Y	0	0
5	IDC ICO H2 -I	N Y	0	0
6	IDC ICO - DO6	N Y	0	0
7	PIN -DO7	N Y	0	0
8	Amptek	N Y	0	0
9		N Y	0	0
10		N Y	0	0
11		N Y	0	0
12		N Y	0	0
13		N Y	0	0
14		N Y	0	0
15		N Y	0	0
16		N Y	0	0

Delay 0.000 (s) Clock 1.000e+007 Hz Update 10.000 Hz

Calcs ENABLE SYNC WITH SCALER+ Less More

5. Depends on experiment, you can open sample stage, IXS/XRS, XES spectrometer from 16idd.adl






## Bluediamond

The Java-based **HPCAT Bluediamond** software has been updated recently, and used together with scanSee program, especially for DCM energy calibration and overlay several scans together.

### Starting Blue Diamond from command line

1. Click Windows icon “start” at the left bottom corner, then click “run”




Every time when this software is started, user should specify the input configuration file “16IDD.txt” in the “C:\HPCATSoftware” directory. Note that this directory is local (not a network directory); however, any arbitrary configuration file in a network directory can be used if properly defined.


The software is straightforward to use and most of the menu items are self-instructing for real-time data display purposes.

## IDL scanSee Software


1. double click “sscan.sav”



2. scanSee window will pop up




3. To start a scan, for example, Sample Z, setup scan range, counting time and number of points, then click "Load&Go"



4. To see the scan, click corresponding detector, for example, here D02 corresponds to IDC ion chamber V1





## 2D scans of Sample Absorption

We want to do a y- and z- scan against a PIN diode to align a sample.

First set up the inner y scan by the scan module and only hit LOAD (NOT Load& Go)

SamY is the 'inside loop' 16IDD:m19.VAL



Counting time (sec) 3,000 1,000 Load&Go Done

Sample Y	Start	End	#Pts	StepSize
Sample Y	-0.200000	0.200000	41	0.010000

Positioner: 16IDD:m19 Units: mm 121

Load RELATIVE AfterScan: PRIOR POS Abort

POSITIONER SETTLING TIME 0.10 (sec)

For IDD scans only !

If you have dummy positioners which are written to the data file, but not scanned  
Fill out the .VAL field with dummies

16IDC:m16.VAL so that the mda file has the correct format to be read by Arun's program

16IDD:scan1 AFTER SCAN: PRIOR POS

I TNFAR	Positioner 1 names	UNITS: mm	Read	Drive	START	CENTER	END	STEP	WIDTH	ARC/PEI
1	16IDC:pm3.RBV	2,500	16IDC:pm3.VAL	0,000	-2,500	0,000	2,500	0,100	5,000	RELATIVE
2	16IDA:m34.RBV	17,357	16IDC:m16.VAL	-0,332	-1,000	0,000	1,000	0,040	2,000	ARSN IITF
3	16IDA:BraggEAO.VAL	10,085	16IDC:m16.VAL	0,000	0,000	0,000	0,000	0,000	0,000	ARSN IITF
4	16IDD:userTran4.K	10,085	16IDC:m16.VAL	0,000	0,000	0,000	0,000	0,000	0,000	ARSN IITF

Make sure the scan is  
RELATIVE

Dummy  
Motors

Then go to the scan2 module

Typically, SamZ is the 'outside loop' 16IDD:m13.VAL

motor 28 is tilt center z for the z in this case  
16IDD:m28.VAL

The screenshot shows the 'scan\_more.adl' window with the following details:

- Title Bar:** scan\_more.adl
- Module:** 2 16IDD:scan2
- Language:** ITIL F
- Feedback:** allback while can record is id
- #PTS:** 41 (highlighted in blue)
- DATA STATE:** POSTFN
- SAVE DATA:** Active
- Scan saved:** 16IDD\_0903.mda
- Positioners:** CLEAR button, SETTLING TIME 0.200 (S)
- Read:** 16IDD:m28.RBV, 4.588
- Drive:** 16IDD:m28.VAL, 3.388
- Parameters:**
  - START: -2.000, CENTER: 0.000, END: 2.000, STEP SIZE: 0.100, WIDTH: 4.000
  - UNITS: mm
  - SCAN MODE: LINEAR
  - ARS/RFI: RELATIVE
  - AFTER SCAN: PRINT PUS
- DetTriggers:** SETTLING TIME 0.200 (S)
- Triggers:** 1 16IDD:scan1.EXSC, 2 (empty)
- Detectors:** 01, 02, 03, 04 (all showing 0.000)
- PLOTS:** Less, More, ? buttons
- Control Buttons:** SCAN (blue), GO (yellow), PAUSE (yellow), ABORT (red)

To execute the scan, hit the SCAN in the scan2 module.

To graph, In scanSee, choose the 'Det 7' and click it, you should see the real time plot

The outside scan position value will show on the top of the plot

## Steering beam procedure

**First tweak** the pitch of the second crystal to optimize the countrate in the ion chamber

The screenshot shows a control window titled "DCMCtl\_hpcat\_new.adl". It contains several sections for beam control:

- DCM beam flux:** IC\_V (130785,00) and IC\_H (115906,50).
- DCM beam steering:** BPM\_V (-2452,40) and BPM\_H (-353,90).
- DCM E (keV):** 7,1320, with a tweak range from -0,0200 to +0,0200 and a "scan" button.
- Pitch tweak section:** A dashed oval highlights the "Pitch" row. It shows a current value of 1098,964500 microns, a target of 1099,000000, and a tweak range from -0,500000 to +0,500000. A "scan" button is next to it.
- Roll tweak section:** Shows a current value of 6137,017875 microns, a target of 6137,000000, and a tweak range from -0,500000 to +0,500000. A "scan" button is next to it.
- Buttons:** "OK", "Done", and a red "All Stop" button.
- DCM message:** A text input field.

scaler16\_more.adl

Done	OneShot	time	Count time	Elapsed time
Count	AutoCount	1.00	1.000	2.895


#	Description	Gate?	Preset count	Actual count
1	Timer	N Y	40000000	28952942
2	IDC IC0 V1 -I	N Y	0	151292
3	IDC IC0 V2 -I	N Y	0	101665
4	IDC IC0 H1 -I	N Y	0	115270
5	IDC IC0 H2 -I	N Y	0	108868
6	IDC IC0 - D0	N Y	0	23167
7	PIN -D07	N Y	0	358
8	AmpTek	N Y	0	3616
9		N Y	0	0
10		N Y	0	0
11		N Y	0	0
12		N Y	0	0
13		N Y	0	0
14		N Y	0	0
15		N Y	0	0
16		N Y	0	0

Delay 0.000 (s) Clock 1.000e+007 Hz Update 10.000 Hz

Calcs ENABLE SYNC WITH SCALER+ Less More

**Find the present vertical and horizontal positions** of the incident beam in IDC

Go to the IDC Slit window




IDC vertical size = 2.5 when taking data

Close size to 0.1 mm

D03 plots the ion chamber counts after the vertical slit

Scan the IDC Vpos (Vertical position) to find the vertical position of the beam




Note center of vertical position scan

then open IDC vertical slit to 2mm, close IDC horizontal to 0.1 (3mm when taking data)

scan horizontal slit position to find the center of horizontal position

if the horizontal position is negative, roll positive



scan horizontal slit position again, adjust roll until position within 50um

then open IDC horizontal slit to 3mm, then close vertical slit to 0.1, scan vertical slit position, if within 50um of position 35, it is ok, otherwise,

.....

Now we will **move the vertical position of the beam** using an IDL program

Open IDL 6.3 icon

Go to 'Files'

Go to 'Recent Files' and open

'DCM\_Vertical\_Vosition.pro'

Compile it (click the icon)

Run it (click the icon)

Input on the command line the vertical position you just found 31.46 return

The macro will move the mono, and scan once again

The position of the beam should move to 35.00 mm (within 50 microns is fine)



If not then iterate by running the IDL program again

When the vertical position is within say 50 microns of 35.000mm:

Make sure the Vsize opens back to 2.000 mm

At this point we have finished moving the beam back to 35.000 at the IDC slits, then **remember to reset the DCM**

.....



Click DCM, select "Full DCM Control", Click "Set", move mouse to energy window "7.112", press "enter", if it shows "done" then ok, otherwise, it shows "moving", then click "move", it will show "Done", then click "OK", then click "Use".




## Steer beam through clean-up Pinhole

Even though the beam is in the correct place in IDC, the focused beam in IDD may not be going thru the cleanup pinhole, so need adjust VKB and HKB pitch by small step, for example, 0.001mrad to get maximum intensity in PID diode (had better remove sample from sample stage)

hpcat08motors.adl

Vertical\_Mirror

VKB upstream moment	3.4650	3.4650	- Tweak +	0.0500	scan	Stop	Move	
	mm					Pause	Go	
VKB downstream moment	4.4751	4.4751	- Tweak +	0.0500	scan	Stop	Move	
	mm					Pause	Go	
Curvature	3.970	3.970	- Tweak +	0.020	scan	Stop	Move	
	mm					Pause	Go	
Ellipticity	-1.010	-1.010	- Tweak +	0.050	scan	Stop	Move	
	mm					Pause	Go	
Height (Sum)	0.161	0.161	- Tweak +	0.020	scan	Stop	Move	
	mm					Pause	Go	
Pitch (Diff)	2.0424	2.0424	- Tweak +	0.0100	scan	Stop	Move	
	mrad					Pause	Go	
Translation (Sum)	-15.000	-15.000	- Tweak +	1.000	scan	Stop	Move	
	mm					Pause	Go	
Rotation (Diff)	0.000	0.000	- Tweak +	0.100	scan	Stop	Move	
	mm					Pause	Go	



**Realign the sample to make sure the beam hit the sample.**


## DCM Energy calibration using Foil

Using Fe foil energy calibration as an example


1. put IDD absorber wheel to -70 degree which is corresponding to Fe foil, if the foil you need is not in the wheel, you can put the wheel to 0 degree then attach your foil to the wheel using tape.




From “Entrance Table” to get Absorber wheel panel, put it at -70 degree for Fe foil, you can click “absorbers” to view correlation between angle and foil.



2. Scan DCM across 7.112keV (Fe k-edge)




Click “DCM” and select “DCM Control New”




Click “scan”, then set up the scan range




Click “Load&Go” to start the DCM energy scan




The red line here shows the DCM energy scan, click “Util” then select “Derivative”



It shows the edge is at 7.1118keV which means the DCM energy is right. If the edge energy is different from tablet edge value, you need reset DCM energy. For example, if the scan shows the edge is at 7.114keV, then you first put DCM energy to 7.114keV



Then select “Full DCM Control” from “DCM” in 16 ID-A control





**kohzuSeqCtl\_All.adl**

H	K	L	E (keV)	lambda (Å)	TH (deg.)	
1	1	1	172.4858	2.35164	34.81790	High Limit
a	3.56688		7.1140	1.74282	25.03366	Actual
2d	4.11868		7.1140	1.74282	25.03366	Desired
			5.2722	0.07188	1.00000	Low Limit
			< 0.0500 >	< 0.010 >	< 0.010 >	

Manual Auto Move

Theta	Y	Z	
16IDA;m34	16IDA;m33	16IDA;m32	Motor PV
25.03366	19.315	41.356	Actual Pos.
25.03366	19.314	41.356	Motor Command
25.03366	19.314	41.356	Ideal
0.05000	0.10000	0.50000	Speed

Channel Cut inhibits Y1 and Z2 motors

Calibration: Use Set y offset: 35.000

Normal Channel Cut Freeze Z Freeze Y

OK

All Stop Done

Click “set” and change 7.114 to 7.112keV

**kohzuSeqCtl\_All.adl**

H	K	L	E (keV)	lambda (Å)	TH (deg.)	
1	1	1	172.4858	2.35164	34.81790	High Limit
a	3.56688		7.1120	1.74331	25.04121	Actual
2d	4.11868		7.1120	1.74331	25.04118	Desired
			5.2722	0.07188	1.00000	Low Limit
			< 0.0500 >	< 0.010 >	< 0.010 >	

Manual Auto Move

Theta	Y	Z	
16IDA;m34	16IDA;m33	16IDA;m32	Motor PV
25.04121	19.316	41.345	Actual Pos.
25.04118	19.316	41.345	Motor Command
25.04118	19.316	41.345	Ideal
0.05000	0.10000	0.50000	Speed

Channel Cut inhibits Y1 and Z2 motors

Calibration: Use Set y offset: 35.000

Normal Channel Cut Freeze Z Freeze Y

OK


All Stop Done



Remember to click “Use” to put DCM into function. You have reset DCM energy, it is better to re-do an edge scan to confirm the edge energy.



## Overplot several scans in Bluediamond


You just finish a scan, for example, the following shows an energy calibration scan. You want to compare it to a similar scan you did before.

1. click “File”, then select “Open”, choose the right scan file to open. Normally the default folder is the folder where the mda file is saved. If not, please select the path “epics on degas/saveData/16idd/run20\*\*\_\*/your folder”,









2. select the file you want to open, for example here, 16IDD\_0070.mda

**HPCAT Old .mda File Selector**

16IDD\_0082.mda
16IDD\_0070.mda


Positioners
Detectors

PV	Description	Minimum	Maximum	Select	Color
16IDD:scaler1.T		1.0	1.0	<input type="radio"/>	...
16IDD:scaler1.S2		61928.0	65332.0	<input type="radio"/>	...
16IDD:scaler1.S3		44615.0	47627.0	<input type="radio"/>	...
16IDD:scaler1.S4		48178.0	52412.0	<input type="radio"/>	...
16IDD:scaler1.S5		55355.0	58358.0	<input type="radio"/>	...
16IDD:scaler1.S6		4649.0	55724.0	<input type="radio"/>	...
16IDD:scaler1.S7		11208.0	24736.0	<input type="radio"/>	...
16IDD:scaler1.S8		5.0	39.0	<input type="radio"/>	...
16IDD:scaler1_calc6.VAL		0.588047...	1.680625...	<input checked="" type="radio"/>	...
16IDD:scaler1_calc8.VAL		-0.69392...	0.682310...	<input type="radio"/>	...
16IDA:m34.REP	DCM Rotation	-186830.0	-182326.0	<input type="radio"/>	...
16IDD:userTran4.M	DCM energy Encoder	2.144712...	2.184747...	<input type="radio"/>	...
S:SRcurrentAI.VAL	SR Current	101.6705...	102.2436...	<input type="radio"/>	...
16IDD:userTran4.K	DCM energy Encoder	7.092016...	7.152036...	<input type="radio"/>	...
16IDA:m34.RCNT	DCM Rotation	0.0	11.0	<input type="radio"/>	...
16IDB:scaler1_cts2.A		27076.66...	27286.66...	<input type="radio"/>	...
16IDD:m18.RBV	Spectrometer Theta	-1.42108...	-1.42108...	<input type="radio"/>	...
16IDD:m27.RBV	Sample Omega	7.644300...	7.644300...	<input type="radio"/>	...

Hide
Show

Close





Here scan82 and scan70 are overplotting.

# Software

## Introduction

Data collection and data analysis software are available to users at ID-D beamline. As a rule please do not install any other software on these computers. If you do need to use your own software for data analysis, please contact the beamline staff and they will work with you to find a solution. Two PCs are provided at each experimental station to assist users in data collection and data analysis. In the next sections, data collection, data processing and general software available on these PCs is listed.

### **Software at Control Computer (16idd):**

1. EPICS WIN32 Extensions
2. Hummingbird Connectivity V8.0
3. IDL 6.3
4. Microsoft Office 2003
5. PowerDVD5.7
6. SSH Secure Shell 3.2.
7. Adobe Reader 7.0
8. Symantec pcAnywhere
9. avast! Antivirus
10. Blue Diamond Version Aug. 17, 2009
11. OriginPro 7.0
12. Roxio (burn CD and DVD)
13. PowderDiffraction (by Wenge Yang)
14. ScanSee
15. MCA

### **Software at Data Analysis Computer (Bella):**

1. EPICS WIN32 Extensions
2. Hummingbird Connectivity V8.0
3. Microsoft Office 2003
4. Roxio
5. IDL 6.2



6. SSH Secure Shell 3.2.9
7. Adobe Reader 8.0
8. Internet Explorer 7.0.5730.11
9. ION 6.4 (unknown distributor)
10. ImageJ 1.42q (seems to be an image processing software)
11. Fit2D V12.077
12. GIMP
13. EditPro (evaluation version)
14. IFeffit (for EXAFS)
15. OriginPro 7.0
16. Acrobat standard 6.0
17. Acrobat distiller 6.0
18. OpenOffice
19. XOP2.3
20. Picasa 2

### **Software at Focus Checking (hpcat36):**

1. Mathematica 4.1
2. Microsoft Office 2003
3. OriginPro 7.0
4. Roxio
5. Sigmaplot 8.0
6. SSH Secure Shell
7. Acrobat standard 6.0
8. Acrobat distiller 6.0
9. WinView32
10. Roper Scientific